

danumed® GastroTube



COMFORT

- Extensive selection of sizes from Fr 12 to Fr 24 for the individual and perfect fit
- Highest patient comfort due to soft, flexible silicone
- Soft and adjustable retention plate for secure ixation
- Structure of the retention plate allows better air circulation in the stoma site area
- Length marking at a distance of 2 cm
- Display of the balloon volume on the shaft
- Side valve access for balloon filling
- Valve colour coding for easy sizing
- Convenient accessories in every single package:
 a pack of lubricant for easy insertion and a danumed[®]
 balloon syringe 10 ml to block the retention balloon.
- Easy to clean and maintain

SAFETY

- Gastric replacement feeding tube made of high-quality medical grade silicone
- All danumed®GastroTubes are free of Latex, PVC, DEHP and BPA
- Flat rounded tip allows atraumatic insertion
- Recessed tip to prevent trauma
- High balloon stability of the inner retention balloon for long usage and secure fixation
- ENFit® connectors according to EN ISO 80369-3 avoid tubing misconnections with other systems
- Available with ENFit® and Y-ENFit® connectors
- Radiopaque stripe along the shaft

danumed [®] GastroTubes with ENFit[®] and Y-ENFit[®]

REF	Product description	Size in Fr	Pieces/ sales unit	Latex- free	PVC- free	DEHP- free	Bisphenol- A-free
DA 1044 0112	danumed®GastroTube with ENFit®	12	1 piece	•	•	•	•
DA 1044 0114	danumed®GastroTube with ENFit®	14	1 piece	•	•	•	•
DA 1044 0116	danumed®GastroTube with ENFit®	16	1 piece	•	•	•	•
DA 1044 0118	danumed®GastroTube with ENFit®	18	1 piece	•	•	•	•
DA 1044 0120	danumed®GastroTube with ENFit®	20	1 piece	•	•	•	•
DA 1044 0122	danumed®GastroTube with ENFit®	22	1 piece	•	•	•	•
DA 1044 0124	danumed®GastroTube with ENFit®	24	1 piece	•	•	•	•
DA 1044 0212	danumed®GastroTube with Y-ENFit®	12	1 piece	•	•	•	•
DA 1044 0214	danumed®GastroTube with Y-ENFit®	14	1 piece	•	•	•	•
DA 1044 0216	danumed®GastroTube with Y-ENFit®	16	1 piece	•	•	•	•
DA 1044 0218	danumed®GastroTube with Y-ENFit®	18	1 piece	•	•	•	•
DA 1044 0220	danumed®GastroTube with Y-ENFit®	20	1 piece	•	•	•	•
DA 1044 0222	danumed®GastroTube with Y-ENFit®	22	1 piece	•	•	•	•
DA 1044 0224	danumed®GastroTube with Y-ENFit®	24	1 piece	•	•	•	•

Each danumed®GastroTube is individually sterile packed. The danumed®GastroTube is delivered in a sterile peel-pouch and additionally with lubricant and a sterile danumed®balloon syringe 10ml in each single box. 10 boxes in a shipping carton.



INNOVATORS IN ENTERAL CARE

For more information, please contact us on (02) 9882 3666 or at amsl@amsl.com.au amsl.com.au

